



Repeatability and Reproducibility of Multiparametric Magnetic Resonance Imaging of the Liver



Henry Wilman^{1,2}, Velicia Bachtiar¹, Jaco Jacobs¹, Rex Newbould¹, Michael Gyngell¹, Cat Kelly¹ Matt Kelly¹, Stefan Neubauer^{1,3} Rajarshi Banerjee¹.

¹Perspectum Diagnostics, Oxford, UK ²University of Westminster, London, UK ³University of Oxford, UK

Background

- ▶ LiverMultiScan™ is a MRI based technology which produces quantitative analysis for hepatic fat, T2* and iron corrected T1 (cT1).
- ▶ cT1 has been shown to correlate with different stages of fibro-inflammatory disease [1] and predicts outcomes [2].
- ▶ LiverMultiScan™ is used in NASH clinical trials and allows the safe and non-invasive characterisation of liver disease.
- ▶ We tested LiverMultiScan™'s metrics for repeatability and reproducibility across different MRI manufacturers, models and software versions and field strengths.

Methodology

- ▶ 61 adult volunteers, a mix of healthy volunteers and those with a background of liver disease. cT1, T2* and Proton Density Fat Fraction maps were acquired on combinations of three Siemens MRI Scanners (1.5T Avanto-Fit; 3T Prisma; 3T Skyra, all VE11C, MyoMaps) and two Phillips MRI scanners (1.5T Ingenia; 3T Ingenia, both 5.3.0, CardiacQuant).
- ▶ Participants were scanned on at least two different scanner models and field strengths in a pseudorandomised order, with a maximum of 1 week between the scans. Two acquisition repeats were conducted for each scan, with participants leaving the scanner between.
- ▶ The same protocol was performed on phantoms across 8 scanner models; the 5 previously mentioned plus 1.5T Phillips Achieva and Phillips AchievaDstream (both 5.3.0, CardiacQuant) and a 1.5T Siemens Aera (VE11C, MyoMaps).
- ▶ Standardisation of T1 maps across the scanner models and software versions was based on 90 previously acquired phantom datasets (figure 1).




Figure 1: Bland-Altman plots demonstrating T1 measures in phantoms (a) before and (b) after standardisation. Reference scanner: 3T Prisma, 1.5T Avanto-fit.




Figure 2: Bland-Altman plots from in vivo measurements across manufacturer and field strength for (a) cT1, (b) T2*, and (c) PDFF. Reference Scanner cT1 and PDFF: 3T Prisma. Reference Scanner T2*: 3T Prisma, 1.5T Avanto-Fit.




Figure 3: Bland-Altman plots from phantom measurements across manufacturer and field strength for (a) T1, (b) T2*, and (c) PDFF. Reference Scanner T1 and PDFF: 3T Prisma. Reference Scanner T1, PDFF, T2*: 3T Prisma, 1.5T Avanto-Fit.

Results

- ▶ Standardised cT1 in participants demonstrated high reproducibility across different scanner models, software versions and field strengths (CoV 3.3%; bias 6.5 ms, 95% LoA of -76.3 ms to 89.2 ms), comparing favourably to MRE's CoV of 10.7% [3] (Figure 2).
- ▶ T2* (CoV 6.6%; bias -1.7ms, 95% LoA of -6.6ms to 3.2 ms) and PDFF measurements (CoV 17%; bias 0.06%, 95% LoA of -0.69% to 0.82%) showed excellent reproducibility across field strengths and scanner models (Figure 2).
- ▶ Bland-Altman analysis of the T1 phantom measurements showed a clear reduction in bias (from -20ms to -4.7ms), tightening of the 95% Limits of Agreement (LoA: from -59.2ms – 19ms, to -25.3ms – 15.9ms) and reduction in mean coefficient of variation (CoV: 2.5% to 1.0%) after standardisation (Figure 3).

Conclusions

- ▶ We demonstrate standardised cT1 is a repeatable and reproducible metric independent of vendor (Philips or Siemens) and field strength (1.5T or 3T).
- ▶ LiverMultiScan™ in full (cT1, T2*, PDFF) represents a robust and reliable non-invasive tool for liver tissue characterisation.

henry.wilman@perspectum-diagnostics.com

- [1] Pavlides M, Banerjee R, Tunnicliffe EM, Kelly C, Collier J, Wang LM, et al. Multi-parametric magnetic resonance imaging for the assessment of non-alcoholic fatty liver disease severity. *Liver International*. 2016; Available from: <http://doi.wiley.com/10.1111/liv.13284>.
- [2] Pavlides M, Banerjee R, Sellwood J, Kelly CJ, Robson MD, Booth JC, et al. Multiparametric magnetic resonance imaging predicts clinical outcomes in patients with chronic liver disease. *Journal of Hepatology*. 2016;64:308–315.
- [3] Trout AT, Serai S, Mahley AD, Wang H, Zhang Y, Zhang B, et al. Liver Stiffness Measurements with MR Elastography: Agreement and Repeatability across Imaging Systems, Field Strengths, and Pulse Sequences. *Radiology*. 2016 December;281(3):793–804.